



HPLC RECORDS

(Ref. SOP LAB-090)

Calibration of HPLC

Date: _____ System: _____ Analyst signature: _____

1. Reproducibility Check

Sample Set Name: _____

Injection Volume: _____ (μL)

Area % RSD: _____ (should be less than 2%)

Retention Time % RSD: _____ (should be less than 1%)

2. Detector Check

(A) Waters 996 Photodiode Array Detector

PDA Internal Diagnostic			
Internal Tests		Pass (P) / Fail (F)	Test display/ Comments
*	CPU		
*	Timer		
*	Shutter		
*	Lamp		
*	Wavelength Accuracy		
*	Communications		
*	Optics		
*	ROM		
*	RAM		
*	DSP		

(B) Waters 470 Scanning Fluorescence Detector

- (i) Check the powering up displays
If: SELF DIAGNOSIS
LAMP EM POOR
is displayed, lamp radiation energy is low and Waters must be called for service.

- (ii) Check the hours of lamp time used.
The lamp has an expected lifetime of 1,000 hours.

LAMP OPERATION TIME _____ HR



HPLC RECORDS

(Ref. SOP LAB-090)

(C) ELSD Detector Check

Diagnostic Test	Current	Max.	Min.	Pass or Fail
1. Optics Test (mV) Laser off: Laser on: Offset:				
2. Nebuliser Gas Pressure Test (psi)				
3. Flow Meter Test (L/min)				

COMMENT: _____

(D) ECD 464 Detector

To ensure good working condition maintenance should be performed every 6 months.
Check chart on ECD detector.

Date working electrode cleaned: _____ Due date: _____

Date changed solution in reference electrode: _____ Due date: _____

Comments: _____

(E) Differential Refractometer 410

Internal Temperature (°C) _____

Auto zero performed _____ (conforms)

Comments: _____



HPLC RECORDS

(Ref. [SOP LAB-090](#))

HPLC Calibration Log

SYSTEM NAME:

Date	HPLC in use- (Yes) or (No)	Date of Calibration	Calibration Pass or Fail	Analyst Signature	Comments



Form-720
Issue date:

HPLC RECORDS

(Ref. [SOP LAB-090](#))

Record of Performance Tests for HPLC Columns

Column Details: Type _____ Length _____ Internal Diameter _____

Supplier's Column Performance Details

Mobile Phase in HPLC Method (M1): _____

Standard Solution used (S1): _____

Product method Details

Product the Column has been dedicated to: _____

Mobile Phase in HPLC Method (M2): _____

Standard Solution used (S2): _____

Date	Analyst	Column Serial Number	HPLC System used	Sample Set name	Standard S1 or S2	Mobile Phase M1 or M2	Flow Rate mL/ min	Injection Volume μ L	λ nm	Plate Counts-N 5 sigma	Capacity Factor k'	USP Tailing