

AUTOCLAVE DRY CYCLE ACCEPTANCE REVIEW CRITERIA

(Ref. SOP VAL-175)

Run Number Cy	ycle Number	Project N	lumber		Autoclave Number	
	A	Oritoria				Acceptance Criteria Achieved (Yes/No/ Accepted Under
Acceptance Criteria Recorded cycle parameters are consistent with programmed cycle parameters.						Deviation)
The Department of Agriculture (DoA) requirements are met for DoA regulated loads (If load core measured, all points ≥121°C for 15min, else all points ≥121°C for 30min and all biological indicators do not show growth).						
The sterilisation phase time is not less than 15 minutes for non DoA regulated loads and temperature regulated cycles.						
The measured chamber pressure is within ± 1°C of the mean chamber temperature at the midpoint of the sterilisation phase. Mean Chamber Temperature =°C Chamber Pressure = kPa Saturation Temperature = °C						
Steam Penetration Studies (Load)						
A minimum accumulated F_0 of 22 minutes is obtained at all load measurement (penetration) locations at the end of the cooling phase on the autoclave printout.						
Min F ₀ Obtained		T	C Location			
All exposed Bls must test negative for growth, except for the positive control.						
The autoclave load is visibly dry. (Note-visibly dry is defined as, upon visual inspection of the load at the end of the cycle, the load items and all layers of wrapping material are found to be free of pooled liquid).						
Heat Distribution Studies (Chamber)						
The measured temperature of the thermocouples within the chamber is within ±1°C of the mean chamber temperature at each time point.						
The measured temperatures within the chamber are within (+2/-1)°C of the						
sterilisation set point.		Se	t Point Temp.			
Lower Temperatur	е	Upper	Temperature			
The autoclave drain probe achieves a minimum temperature of ≥1°C below the						
sterilisation set point.	Set Point		Min. T	emp.		
Test Acceptance	☐ Pass ☐ Fail					
Comments						
	1					
Tested by/Date			Reviewed b	y/Date		

Form 835 Issue date: Version: 1 Page: 1 of 1