

Department	Micro Laboratory	Document no	MICLAB – METHOD 017		
Title	Method for Preparation of DNase Agar (Oxoid CM321)				
Prepared by:		Date:		Supersedes:	
Checked by:		Date:		Date Issued:	
Approved by:		Date:		Review Date:	

1.0 **SUMMARY OF CHANGES**

Version #	Revision History
MICLAB – METHOD 017	New

2.0 **PURPOSE**

This document describes the method for preparation of DNase Agar (Oxoid CM321).

3.0 **SCOPE**

DNase Agar is used for the detection of microbial deoxyribonuclease enzymes, particularly from Staphylococci. The DNase reaction for Staphylococci is an indication of pathogenicity; it cannot be used as the sole criteria for identification.

4.0 **RESPONSIBILITY \ BUSINESS RULES**

All microbiology staff at GMP site.

5.0 **PROCEDURE**

5.1 Materials and Reagents required

- 5.1.1 Plastic spoon
- 5.1.2 Measuring Cylinder
- 5.1.3 RO Water
- 5.1.4 Microwaveable beaker
- 5.1.5 Plastic Pipettes
- 5.1.6 DNase Agar powder (Oxoid CM321)
- 5.1.7 Clean Glassware (flat medical or 250mL schott bottles)
- 5.1.8 Sterile Petri dishes
- 5.1.9 1N NaOH and 1N HCl

5.2 Method

- 5.2.1 Weigh out 39g per litre of DNase Agar powder into a beaker.
- 5.2.2 Add required volume of RO water and mix well.

Department	Micro Laboratory	Document no	MICLAB – METHOD 016
Title	Method for Preparation of DNase Agar (Oxoid CM321)		

- 5.2.3 Dissolve on high setting in microwave with frequent stirring.
- 5.2.4 Adjust pH to 7.5 using 1N NaOH or 1N HCl.
- 5.2.5 Dispense 150mL or 300mL volumes into flat medical bottles or schott bottles respectively.
- 5.2.6 Sterilise by autoclaving at 121°C for 15 minutes.
- 5.2.7 Final pH after sterilisation should be 7.3 ± 0.2 .
- 5.2.8 Mix gently before pouring into sterile petri dishes.
- 5.2.9 Record all details of media preparation on SF150712.

5.3 Quality Control Requirements

QUALITY CONTROL REQUIREMENTS			
STORAGE:	BULK	-	6 months in dark cupboard
	POURED PLATES	-	2 weeks in fridge (2-8°C)
ECOMETRIC EVALUATION / FERTILITY			
CONTROL ORGANISMS:			
Positive:	Staphylococcus aureus N.C.T.C 6571	Negative:	N/A
Growth Index:	≥ 3	Growth Index:	N/A
INCUBATION CONDITIONS:			
Temperature:	37 ± 1°C	Time:	18 hours

6.0 DEFINITIONS / ACRONYMS

NA