

<b>Department</b>	<b>Micro Laboratory</b>	<b>Document no</b>	<b>MICLAB – METHOD 031</b>	
<b>Title</b>	<b>Method of preparation of Bacteriological Peptone (Oxide L34)</b>			
Prepared by:		Date:		Supersedes:
Checked by:		Date:		Date Issued:
Approved by:		Date:		Review Date:

## 1.0 **SUMMARY OF CHANGES**

<b>Version #</b>	<b>Revision History</b>
<b>MICLAB METHOD 031</b> –	New

## 2.0 **PURPOSE**

This document describes the method of preparation of Bacteriological Peptone (Oxide L34).

## 3.0 **SCOPE**

Bacteriological peptone (0.1%) is used as a high quality diluent for microbiology testing.

## 4.0 **RESPONSIBILITY \ BUSINESS RULES**

All microbiology staff at Pfizer, Caringbah.

## 5.0 **PROCEDURE**

### **5.1 Materials and Reagents Required**

- 5.1.1 Bacteriological Peptone Neutralised (Oxoid L34)
- 5.1.2 Beaker
- 5.1.3 Measuring cylinder
- 5.1.4 1N NaOH and 1N HCl
- 5.1.5 Plastic pipettes
- 5.1.6 RO water
- 5.1.7 Automatic dispenser or pouring jug and funnel
- 5.1.8 Plastic spoon

### **5.2 Method**

- 5.2.1 Add 1g per litre of Neutralised Bacteriological Peptone into approximately 500mL of RO water, mix to dissolve. Make up to 1 litre with RO water.
- 5.2.2 Adjust pH to 7.0 and dispense into required glassware (eg. 9mL & 9.9mL volumes). If precise volume is required use automatic dispenser.
- 5.2.3 Sterilise by autoclaving at 121°C for 15 minutes.
- 5.2.4 pH after autoclaving 7.0 ± 0.2.
- 5.2.5 Record all details of media preparation on SF150712.

### **5.3 Quality Control Requirements**

<b>Department</b>	<b>Micro Laboratory</b>	<b>Document no</b>	<b>MICLAB – METHOD 031</b>
<b>Title</b>	<b>Method for preparation of Pseudomonas Agar Base (Oxoid CM 559)</b>		

<b>QUALITY CONTROL REQUIREMENTS</b>			
STORAGE:	BULK	-	6 months in dark cupboard
<b>ECOMETRIC EVALUATION / FERTILITY</b>			
CONTROL ORGANISMS:			
Positive:	Pseudomonas aeruginosa	Negative:	Staphylococcus aureus
Growth Index:	Growth	Growth Index:	Growth
INCUBATION CONDITIONS:			
Temperature:	30 ± 1°C	Time:	24 hours

For sampling Plan and Acceptance Criteria for microbiological culture media refer to SF150154.

## **6.0 DEFINITIONS / ACRONYMS**

N/A.

## **7.0 REFERENCES**

7.1 The Oxoid Manual 8<sup>th</sup> Edition 1998.